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1. Reichel, H., Koeffler, H.P., Norman, A.W. The role of the vitamin D endocrine system in health and disease. N. Engl. J. Med. 320: 981-991. (1989)
2. Holick, M.F., Vitamin D: Photobiology, Metabolism, etc., In "Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism", American Society for Bone and Mineral Research., 3rd Ed. Lippincott-Raven, Philadelphia, 74-81. (1996).
3. Chapuy M-C, Preziosi P, Maamer M, et al. Prevalence of vitamin D insufficiency in an adult normal population. Osteoporos Int 1997;7: 439-443.
4. Heaney, RP. Functional indices of vitamin D status and ramifications of vitamin D deficiency. Am J Clin Nutr 2004;80 (suppl):1706S-1709S.
5. Holick MF. Sunlight and vitamin D for bone health and prevention of autoimmune diseases, cancers, and cardiovascular disease. Am J Clin Nutr 2004; 80: 1678S - 1688S.
6. Webb AR, Kline L, Holick MF. Influence of season and latitude on the cutaneous synthesis of vitamin D3: exposure to winter sunlight in Boston and Edmonton will not promote vitamin D3 synthesis in human skin. J Clin Endocrinol Metab 1988;67: 373-378.
7. Holick MF, Matsuoka LY, Wortsman J. Age, vitamin D and solar ultraviolet. Lancet 1989;ii:1104-1105.
8. Matsuoka LY, Iide L, Wortsman J et al. Sunscreens suppress cutaneous vitamin D3 synthesis. J Clin Endocrinol Metab 1987;64:1165-1168.
9. Salih FM. Effect of clothing varieties on solar photosynthesis of previtamin D3: an in vitro study. Photodermatol Photoimmunol Photomed 2004;20: 53-58.
10. Clemens TL, Henderson SL, Adams JS, Holick MF. Increased skin pigment reduces the capacity of skin to synthesize vitamin D3. Lancet 1982;9:74-76.
11. National Osteoporosis Foundation. Prevention -Vitamin D. <http://www.nof.org/prevention/vitaminD.htm>
12. KDOQI Clinical Practice Guidelines for Bone Metabolism and Disease in Children With Chronic Kidney Disease. [http://www.kidney.org/PROFESSIONALS/kdoqi/guidelines\\_pedbone/guide8.htm](http://www.kidney.org/PROFESSIONALS/kdoqi/guidelines_pedbone/guide8.htm) Doc: AC-57PL-A  
Issue: 8  
2009年8月4日

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### ab181421 Human TNF alpha SimpleStep ELISA® Kit

For the quantitative measurement of TNF alpha in human serum, plasma, and cell culture supernatant.  
For research use only - not intended for diagnostic use.

For overview, typical data and additional information please visit: [www.abcam.com/ab181421](http://www.abcam.com/ab181421)

**Storage and Stability:** Store kit at +4°C immediately upon receipt. Refer to list of materials supplied for storage conditions of individual components. Observe the storage conditions for individual prepared components in the Materials Supplied section.

#### Materials Supplied

Item	Quantity	Storage Condition
Human TNF alpha Capture Antibody 10X	600 $\mu$ L	+4°C
Human TNF alpha Detector Antibody 10X	600 $\mu$ L	+4°C
Human TNF alpha Lyophilized Recombinant Protein	2 Vials	+4°C
Antibody Diluent 4BR	6 mL	+4°C
10X Wash Buffer PT	20 mL	+4°C
TMB Development Solution	12 mL	+4°C
Stop Solution	12 mL	+4°C
Sample Diluent NS	50 mL	+4°C
SimpleStep Pre-Coated 96-Well Microplate	96 Wells	+4°C
Plate Seal	1	+4°C

#### Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully utilize this assay:  
Microplate reader capable of measuring absorbance at 450 or 600 nm.  
Method for determining protein concentration (BCA assay recommended).  
Deionized water.  
PBS [1.4 mM KH<sub>2</sub>PO<sub>4</sub>, 8 mM Na<sub>2</sub>HPO<sub>4</sub>, 140 mM NaCl, 2.7 mM KCl, pH 7.4].  
Multi- and single-channel pipettes.  
Tubes for standard dilution.  
Plate shaker for all incubation steps.  
Optional: Phenylmethylsulfonyl Fluoride (PMSF) (or other protease inhibitors).

#### Reagent Preparation

Equilibrate all reagents to room temperature (18-25°C) prior to use. The kit contains enough reagents for 96 wells. **The sample volumes below are sufficient for 48 wells (6 x 8-well strips); adjust volumes as needed for the number of strips in your experiment.** Prepare only as much reagent as is needed on the day of the experiment. Capture and Detector Antibodies have only been tested for stability in the provided 10X formulations.

#### Sample Preparation

TYPICAL SAMPLE DYNAMIC RANGE	
Sample Type	Range
Serum*	$\leq$ 1:40
Plasma – Citrate*	$\leq$ 1:40
Plasma – EDTA*	$\leq$ 1:40
Plasma – Herapin*	$\leq$ 1:40
PBMC cell culture supernatant**	1:80 -1:10

\*Based on spiked sample  
\*\*Range can vary depending on stimulation factors

**Plasma** Collect plasma using citrate, EDTA or heparin. Centrifuge samples at 2,000 x g for 10 minutes. Dilute samples at least 1:40 into Sample Diluent NS and assay. Store undiluted plasma samples at -20°C or below for up to 3 months. Avoid repeated freeze thaw cycles.

**Serum** Samples should be collected into a serum separator tube. After clot formation, centrifuge samples at 2,000 x g for 10 minutes and collect serum. Dilute samples at least 1:40 into Sample Diluent NS and assay. Store un-diluted serum at -20°C or below. Avoid repeated freeze-thaw cycles.

**Cell Culture Supernatants** Centrifuge cell culture media at 2,000 x g for 10 minutes to remove debris. Collect supernatants and assay. Or dilute samples at least 1:10 into Sample Diluent NS and assay. Store un-diluted samples at -20°C or below. Avoid repeated freeze thaw cycles.

#### Plate Preparation

- The 96 well plate strips included with this kit are supplied ready to use. It is not necessary to rinse the plate prior to adding reagents.
- Unused plate strips should be immediately returned to the foil pouch containing the desiccant pack, resealed and stored at 4°C.
- For each assay performed, a minimum of two wells must be used as the zero control. For statistical reasons, we recommend each sample should be assayed with a minimum of two replicates (duplicates).
- Differences in well absorbance or "edge effects" have not been observed with this assay.

#### Assay Procedure

Equilibrate all materials and prepared reagents to room temperature prior to use. We recommend that you assay all standards, controls and samples in duplicate.

- Prepare all reagents, working standards, and samples as directed in the previous sections.

**1X Wash Buffer PT:** Prepare 1X Wash Buffer PT by diluting 10X Wash Buffer PT with deionized water. To make 50 mL 1X Wash Buffer PT combine 5 mL 10X Wash Buffer PT with 45 mL deionized water. Mix thoroughly and gently.

**Antibody Cocktail:** Prepare Antibody Cocktail by diluting the capture and detector antibodies in Antibody Diluent 4BR. To make 3 mL of the Antibody Cocktail combine 300  $\mu$ L 10X Capture Antibody and 300  $\mu$ L 10X Detector Antibody with 2.4 mL Antibody Diluent 4BR. Mix thoroughly and gently.

#### Standard Preparation

Prepare serially diluted standards immediately prior to use. Always prepare a fresh set of positive controls for every use. The following section describes the preparation of a standard curve for duplicate measurements (recommended).

- IMPORTANT:** If the protein standard vial has a volume identified on the label, reconstitute the TNF alpha standard by adding that volume of Sample Diluent NS indicated on the label. Alternatively, if the vial has a mass identified, reconstitute the TNF alpha standard by adding 1000  $\mu$ L Sample Diluent NS. Hold at room temperature for 10 minutes and mix gently. This is the 10000 pg/mL Stock Standard Solution.
- Label eight tubes, Standards 1–8.
- Add 360  $\mu$ L Sample Diluent NS into tube number 1 and 150  $\mu$ L of Sample Diluent NS into numbers 2–8.
- Use the Stock Standard to prepare the following dilution series. Standard #8 contains no protein and is the Blank control:

Standard #	Dilution Sample	Volume to Dilute ( $\mu$ L)	Volume of Diluent ( $\mu$ L)	Starting Conc. (pg/mL)	Final Conc. (pg/mL)
1	Stock	40	360	10,000	1,000
2	Standard#1	150	150	1,000	500
3	Standard#2	150	150	500	250
4	Standard#3	150	150	250	125
5	Standard#4	150	150	125	62.5
6	Standard#5	150	150	62.5	31.25
7	Standard#6	150	150	31.25	15.63
8	Blank Control	0	300	0	0

To convert sample values obtained with the kit to approximate NIBSC 12/154 units, use the following equation: NIBSC (12/154) approximate value (IU/mL) = 0.094 x SimpleStep Human TNF alpha value (pg/mL).

- Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, resealed and return to 4°C storage.
- Add 50  $\mu$ L of all sample or standard to appropriate wells.
- Add 50  $\mu$ L of the Antibody Cocktail to each well.
- Seal the plate and incubate for 1 hour at room temperature on a plate shaker set to 400 rpm.
- Wash each well with 3 x 350  $\mu$ L 1X Wash Buffer PT. Wash by aspirating or decanting from wells then dispensing 350  $\mu$ L 1X Wash Buffer PT into each well. Complete removal of liquid at each step is essential for good performance. After the last wash invert the plate and blot it against clean paper towels to remove excess liquid.
- Add 100  $\mu$ L of TMB Development Solution to each well and incubate for 10 minutes in the dark on a plate shaker set to 400 rpm. Given variability in laboratory environmental conditions, optimal incubation time may vary between 5 and 20 minutes. Note: The addition of Stop Solution will change the color from blue to yellow and enhance the signal intensity about 3X. To avoid signal saturation, proceed to the next step before the high concentration of the standard reaches a blue color of O.D.600 equal to 1.0.
- Add 100  $\mu$ L of Stop Solution to each well. Shake plate on a plate shaker for 1 minute to mix. Record the OD at 450 nm. This is an endpoint reading.
- Alternative to 7–8: Instead of the endpoint reading at 450 nm, record the development of TMB Substrate kinetically. Immediately after addition of TMB Development Solution begin recording the blue color development with elapsed time in the microplate reader prepared with the following settings:

Mode	Kinetic
Wavelength	600 nm
Time	up to 20 min
Interval	20 sec - 1 min
Shake	Shake between readings

**Note** that an endpoint reading can also be recorded at the completion of the kinetic read by adding 100  $\mu$ L Stop Solution to each well and recording the OD at 450 nm.

#### Download our ELISA guide for technical hints, results, calculation, and troubleshooting tips:

[www.abcam.com/protocols/the-complete-elisa-guide](http://www.abcam.com/protocols/the-complete-elisa-guide)

For technical support contact information, visit: [www.abcam.com/contactus](http://www.abcam.com/contactus)

Copyright © 2021 Abcam. All Rights Reserved. All information / detail is correct at time of going to print.  
Version 4 | 2021-11-08

## ab181421 Human TNF alpha SimpleStep ELISA® Kit

### ASSAY SPECIFICITY

This kit is designed for the quantification of Human TNF alpha.

Native signal we detected in cell culture supernatant

Spiked protein experiments were used to validate serum, plasma citrate, plasma EDTA, plasma heparin and cerebrospinal fluid

### CROSS REACTIVITY –

Recombinant proteins were prepared at 1000 pg/mL and assayed for cross reactivity. No cross reactivity was found for the following targets: - Human IL-2 - Human IL-4 - Human IL-1 alpha - Human IFN gamma - Human TNF beta - Human TNF R1.

### INTERFERENCE –

Recombinant Human TNF R1 was prepared at 1000 pg/mL and tested for interference. No interference with was observed.

### SPECIES REACTIVITY -

This kit recognizes human TNF alpha protein.

Other species reactivity was determined by measuring 1000 pg/mL recombinant proteins of various species, interpolating the protein concentrations from the human standard curve, and expressing the interpolated concentrations as a percentage of the protein concentration assayed at the same dilution. Reactive species: Primate

Reactivity < 3% was determined for the following species: Mouse/Rat

### CALIBRATION

This immunoassay is calibrated against a highly purified human TNF alpha. The NIBSC/WHO unclassified purified human TNF alpha preparation 12/154 was evaluated in this kit. The dose response curve of the unclassified standard TNF alpha parallels the SimpleStep standard curve. To convert sample values obtained with the SimpleStep Human TNF alpha kit to approximate NIBSC 12/154 International units, use the equation below.

NIBSC (12/154) approximate value (IU/mL) = 0.094 x SimpleStep Human TNF alpha value (pg/mL).

### CALCULATION –

- Calculate the average absorbance value for the blank control (zero) standards. Subtract the average blank control standard absorbance value from all other absorbance values.
  - Create a standard curve by plotting the average blank control subtracted absorbance value for each standard concentration (y-axis) against the target protein concentration (x-axis) of the standard. Use graphing software to draw the best smooth curve through these points to construct the standard curve.
- Δ Note: Most microplate reader software or graphing software will plot these values and fit a curve to the data. A four parameter curve fit (4PL) is often the best choice; however, other algorithms (e.g. linear, semi-log, log/log, 4 parameter logistic) can also be tested to determine if it provides a better curve fit to the standard values.
- Determine the concentration of the target protein in the sample by interpolating the blank control subtracted absorbance values against the standard curve. Multiply the

resulting value by the appropriate sample dilution factor, if used, to obtain the concentration of target protein in the sample.

- Samples generating absorbance values greater than that of the highest standard should be further diluted and reanalyzed. Similarly, samples which measure at an absorbance values less than that of the lowest standard should be retested in a less dilute form.

### TYPICAL DATA

Typical standard curve – data provided for demonstration purposes only. A new standard curve must be generated for each assay performed

Concentration (pg/mL)	Standard Curve Measurement		Mean O.D
	1	2	
0	0.061	0.068	0.065
15.63	0.127	0.123	0.125
31.25	0.189	0.206	0.197
62.5	0.312	0.312	0.312
125	0.555	0.521	0.538
250	1.065	1.030	1.048
500	1.925	1.941	1.933
1000	3.362	3.161	3.262

Figure 1. Example of human TNF alpha standard curve in Sample Diluent NS. The TNF alpha standard curve was prepared as described in Section 10. Raw data values are shown in the table. Background-subtracted data values (mean +/- SD) are graphed.

### TYPICAL SAMPLE VALUES

#### Sensitivity:

The calculated minimal detectable dose (MDD) is 4.32 pg/mL. The MDD was determined by calculating the mean of zero standard replicates (n=24) and adding 2 standard deviations then extrapolating the corresponding concentration.

### Recovery

Three concentrations of TNF alpha recombinant protein were spiked in duplicate to the indicated biological matrix to evaluate signal recovery in the working range of the assay.

Sample Type	Average Recovery	% Range (%)
2.5% Serum	88	83-93
2.5% Plasma-Citrate	91	86-97
2.5% Plasma-EDTA	100	95-104
2.5% Plasma-Heparin	98	89-110
5% PHA-M treated PBMC Cell Culture supernatant*	110	105-113
50% Cerebrospinal Fluid	83	80 - 88

\*Media is RPMI 1640 containing 10% fetal calf serum.

### Linearity of Dilution

Linearity of dilution is determined based on interpolated values from the standard curve. Linearity of dilution defines a sample concentration interval in which interpolated target concentrations are directly proportional to sample dilution.

Native TNF alpha was measured in the following biological samples in a 2-fold dilution series. Sample dilutions are made in Sample Diluent NS.

Dilution Factor	Interpolated value	10% PHA-M Treated PBMC Cell Culture Supernatant
Undiluted	pg/mL	726.56
	% Expected value	100
2	pg/mL	355.55
	% Expected value	98
4	pg/mL	183.66
	% Expected value	101
8	pg/mL	94.78
	% Expected value	104

Recombinant TNF alpha was spiked into the following biological samples and diluted in a 2-fold dilution series in Sample Diluent NS.

Dilution Factor	Interpolated value	2.5% Human Serum	2.5% Human Plasma (Citrate)	2.5% Human Plasma (EDTA)	2.5% Human Plasma (Heparin)	50% Human Cerebro-spinal Fluid
Undiluted	pg/mL	447.20	451.09	483.71	499.11	392.65
	% Expected value	100	100	100	100	100
2	pg/mL	223.35	234.34	244.59	237.17	192.55
	% Expected value	100	104	101	95	98
4	pg/mL	119.50	127.27	123.94	120.73	110.67
	% Expected value	107	113	102	97	113
8	pg/mL	60.23	59.54	62.05	60.11	56.51
	% Expected value	108	106	103	96	115
16	pg/mL	28.78	NL	NL	NL	NL
	% Expected value	103	NL	NL	NL	NL

Twenty individual healthy human female/male donors were measured in duplicate for the presence of TNF alpha. All values were below the detectable range of the assay NL – Non-Linear

### Precision

Mean coefficient of variations of interpolated values of TNF alpha from a single concentration of PHA-M treated PBMC cell culture supernatant within the working range of the assay.

	Intra-assay	Inter-assay
N=	8	3
CV (%)	2.5	3.1

Download our ELISA guide for technical hints, results, calculation, and troubleshooting tips:

[www.abccom.com/protocols/the-complete-elisa-guide](http://www.abccom.com/protocols/the-complete-elisa-guide)

For technical support contact information, visit: [www.abccom.com/contactus](http://www.abccom.com/contactus)

Copyright © 2021 Abcam. All Rights Reserved. All information / detail is correct at time of going to print. Version 4a | 2021-11-08

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### ab287806 – Kanamycin ELISA Kit

For the quantitative measurement of Kanamycin in Tissue, Milk and milk powder.  
For research use only - not intended for diagnostic use.

For overview, typical data and additional information please visit:  
<http://www.abcam.com/ab287806>

#### Storage and Stability

The entire kit may be stored at 4°C for up to 12 months from the date of shipment. Opened kit may be stable for 1 month at 4°C.

#### Materials Supplied

Item	Quantity	Storage Condition
Micro ELISA Plate	8 x 12 Strips	4°C
Standard (S0-S5)	1 mL x 6	4°C
HRP-Conjugate	7 mL	4°C
Antibody	7 mL	4°C
TMB Substrate	12 mL	4°C
Stop Solution	10 mL	4°C
Sample Diluent	20 mL	4°C
10X Wash Buffer	30 mL	4°C
Plate sealers	4 units	4°C

#### Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully utilize this assay:

- Chemicals: Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O, NaH<sub>2</sub>PO<sub>4</sub>·2H<sub>2</sub>O, and Trichloroacetic acid (TCA)
- Microplate reader capable of measuring absorbance at 450 nm
- 25°C incubator
- Precision pipettes with disposable tips
- Distilled or deionized water
- Clean eppendorf tubes and graduated cylinders for preparing standards or sample dilutions
- Absorbent paper

#### Reagent Preparation

- Bring all reagents to room temperature (20-25°C) 30 minutes before use
- Before using the kit, spin tubes and bring down all components to the bottom of tubes.

**Extraction Solution 1:** Weigh 5.37 g of Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O and 0.78 g of NaH<sub>2</sub>PO<sub>4</sub>·2H<sub>2</sub>O to 100 ml of deionized water, mix thoroughly.

**Extraction Solution 2:** Dilute 5 ml of Sample Diluent into 95 ml deionized or distilled water, mix well.

**Extraction Solution 3:** Dilute 3 g of TCA into 100 ml deionized or distilled water, shake well

- or autowasher, and let it stand for 30 seconds, complete removal of liquid at each step is essential to good performance.
4. Add 100 µl of TMB Substrate to each well, mix well, incubate for 15 minutes at 25°C. Protect from light.
  5. Add 50 µl of Stop Solution to each well, gently tap the plate to ensure thorough mixing.
  6. Read result at 450 nm within 10 minutes

#### Calculations

The mean values of the absorbance values obtained for the standards and the samples are divided by the absorbance value of the first standard (zero standard) and multiplied by 100%. The zero standard is thus made equal to 100% and the absorbance values are quoted in percentages.

$$\text{Absorbance Value (\%)} = B/B0 \times 100\%$$

B: The average absorbance value of the sample or standard  
B<sub>0</sub>: The average absorbance value of the 0 ppb standard

To draw a standard curve: Take the absorbency value of standards as y-axis, logarithmic of the concentration of the Kanamycin standards solution (ppb) as x-axis. The Kanamycin concentration of each sample (ppb), which can be read from the calibration curve, is multiplied by the corresponding dilution factor of each sample, allowing for the actual concentration of sample to be obtained.

Download our ELISA guide for technical hints, results, calculation, and troubleshooting tips:  
[www.abcam.com/protocols/the-complete-elisa-guide](http://www.abcam.com/protocols/the-complete-elisa-guide)

For technical support contact information, visit: [www.abcam.com/contactus](http://www.abcam.com/contactus)

Copyright © 2021 Abcam. All Rights Reserved. All information / detail is correct at time of going to print. Version 1 | 2021-September-14

**Wash Buffer (1X):** If crystals have formed in the concentrate, warm up to room temperature and mix gently until the crystals have completely dissolved. Dilute 10 mL of 10X Wash Buffer into 90 mL deionized or distilled water to prepare 100 ml of 1X Wash Buffer. Keep it at 4°C for one month.

#### Standard Preparation

Standards	S0	S1	S2	S3	S4	S5
Concentration (ppb)	0	0.5	1.5	4.5	13.5	40.5

#### Sample Preparation

- The prepared sample maybe stored for up to one day at 2-8°C

#### Milk and Fresh Milk:

1. Bring the milk sample to room temperature.
2. Transfer 100 µL of sample into a new centrifugal tube and add 900 µL of Extraction Solution 1, shake well.
3. Take 50 µL of sample for further analysis.  
**Δ Note:** Dilution factor: 10

#### Milk Powder:

1. Weigh 1 g of milk powder sample, add 5ml of Extraction Solution 2, shake properly for 5 min.
2. Add 4 ml of Extraction Solution 3, shake properly for 5 min. Centrifuge at 4000 rpm for 10 min.
3. Transfer 100 µL of supernatant into a new centrifugal tube, add 200 µL of Sample Diluent, shake well.
4. Take 50 µL of sample for further analysis.  
**Δ Note:** Dilution factor: 30

#### Tissue (chicken, pork):

1. Weigh 1 g of the homogenized tissue sample, add 10 ml of Extraction Solution 1, shake properly for 5 min. Centrifuge at 4000 rpm for 10 min.
2. Take 50 µL of supernatant sample for further analysis  
**Δ Note:** Dilution factor: 10

#### Assay Procedure

- Bring all reagents and samples to room temperature 30 minutes prior to the assay.
- It is recommended that all standards and samples be run at least in duplicate.
- A standard curve should be run for each assay.

1. Prepare all reagents, samples, and standards.
2. Add 50 µL of Standard or Sample per well. Then add 50 µL of HRP-conjugate to each well and 50 µL of Antibody to each well. Cover the microtiter plate with a new adhesive strip and mix well, then incubate for 30 min at 25°C.
3. Aspirate each well and wash, repeating the process 4 times. Wash by filling each well with 250 µL of Wash Buffer using a squirt bottle, multi-channel pipette, manifold dispenser,

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Version 2a Last updated 10 May 2021

# ab274406 – Bovine Lactoferrin ELISA Kit

For the determination of Lactoferrin in bovine biological samples.

View ab274406

Bovine Lactoferrin ELISA Kit datasheet:

[www.abcam.com/ab274406](http://www.abcam.com/ab274406)

(use [www.abcam.cn/ab274406](http://www.abcam.cn/ab274406) for China, or [www.abcam.co.jp/ab274406](http://www.abcam.co.jp/ab274406) for Japan)

This product is for research use only and is not intended for diagnostic use.

Copyright © 2021 Abcam. All rights reserved


## Table of Contents

1. Overview	1
2. Protocol Summary	2
3. Precautions	3
4. Storage and Stability	3
5. Limitations	4
6. Materials Supplied	4
7. Materials Required, Not Supplied	5
8. Technical Hints	6
9. Reagent Preparation	7
10. Sample Collection and Storage	8
11. Sample Preparation	9
12. Standard Curve preparation	10
13. Assay Procedure	11
14. Calculations	12
15. Notes	13
<b>Technical Support</b>	<b>14</b>

## 1. Overview

The Bovine Lactoferrin ELISA kit (ab274406) is a highly sensitive two-site enzyme linked immunoassay (ELISA) for measuring Lactoferrin (LF) in bovine biological samples.

## 2. Protocol Summary

### **3. Precautions**

**Please read these instructions carefully prior to beginning the assay.**

All kit components have been formulated and quality control tested to function successfully as a kit. Modifications to the kit components or procedures may result in loss of performance.

### **4. Storage and Stability**

- Store the kit at +4 to +8°C in the dark immediately upon receipt.
- Kit has a storage time of 1 year from receipt, providing components have not been reconstituted.
- Refer to list of materials supplied for storage conditions of individual components. Observe the storage conditions for individual prepared components in the Reagent Preparation section.

## 5. Limitations

- This assay will perform as described only when the assay procedure is carefully followed and with adherence to good laboratory practice.
- Factors that might affect the performance of the assay include proper instrument function, cleanliness of glassware, quality of distilled or deionized water, accuracy of reagent and sample pipetting, washing technique, incubation time and/or temperature.
- Do not mix or substitute reagents with those from other lots or sources.

## 6. Materials Supplied

Item	Quantity	Storage Condition
Antibody coated plate	96 wells	4°C
Bovine Lactoferrin Calibrator (Lyophilized)	1 vial	4°C
5X Diluent Buffer concentrate	50 mL	4°C
20X Wash Buffer Concentrate	50 mL	4°C
100X HRP Conjugated Bovine Lactoferrin antibody	1 vial	4°C (Store in dark)
Chromogen Substrate Solution	12 mL	4°C (Store in dark)
Stop Solution	12 mL	4°C

## 7. Materials Required, Not Supplied

These materials are not included in the kit, but will be required to successfully perform this assay:

- Microplate reader capable of measuring absorbance at 450 nm.
- Microplate washer/ aspirator
- Precision pipettes to deliver 2  $\mu$ L to 200  $\mu$ L volumes.
- Adjustable 1-25 mL pipettes for reagent preparation.
- Timer
- Distilled or deionized water.
- Log-log graph paper or computer and software for ELISA data analysis.
- Tubes to prepare standard or sample dilutions.
- Assorted glassware for the preparation of reagents and buffer solutions.

## 8. Technical Hints

- Samples which generate values that are greater than the most concentrated standard should be further diluted in the appropriate sample dilution buffer.
- Avoid foaming or bubbles when mixing or reconstituting components.
- Avoid cross contamination of samples or reagents by changing tips between sample, standard and reagent additions.
- Ensure plates are properly sealed or covered during incubation steps.
- Completely aspirate all solutions and buffers during wash steps. When preparing your standards, it is critical to briefly spin down the vial first. The powder may adhere to the cape and not be included in the standard solution resulting in an incorrect concentration. Be sure to dissolve the powder thoroughly when reconstituting. After adding Assay Diluent to the vial, we recommend inverting the tube a few times, then flick the tube a few times, and then spin it down; repeat this procedure 3-4 times. This is an effective technique for thorough mixing of the standard without using excessive mechanical force.
- Do not vortex the standard during reconstitution, as this will destabilize the protein.
- Once your standard has been reconstituted, it should be used right away or else frozen for later use.
- Keep the standard dilutions on ice during preparation, but the ELISA procedure should be done at room temperature.
- Be sure to discard the working standard dilutions after use – they do not store well.
- This kit is sold based on number of tests. A 'test' simply refers to a single assay well. The number of wells that contain sample, control or standard will vary by product. Review the protocol completely to confirm this kit meets your requirements. Please contact our Scientific Support staff with any questions.

## 9. Reagent Preparation

- Equilibrate all reagents to room temperature (18-25°C) prior to use. The kit contains enough reagents for 96 wells.

### 9.1 1X Diluent Solution

The diluent solution is supplied as 5X Diluent Concentrate and must be diluted 1/5 with distilled or deionized water (1 part buffer concentrate, 4 parts dH<sub>2</sub>O). The 1X Diluent Solution is stable for at least one week from the date of preparation and should be stored at +4-8°C.

### 9.2 1X Wash Buffer

The wash solution is supplied as 20X Concentrate and must be diluted 1/20 with distilled or deionized water (1 part buffer concentrate, 19 parts dH<sub>2</sub>O). Crystal formation in the concentrate is not uncommon when storage temperatures are low. Warming of the concentrate to 30 - 35°C before dilution can dissolve crystals. The 1X Wash Buffer is stable for at least one week from the date of preparation and can be stored at +4-8°C.

### 9.3 1X HRP-Conjugated Bovine lactoferrin antibody

Calculate the required amount of 1X Enzyme-Antibody Conjugate solution for each microtitre plate test strip by adding 10 µL Enzyme-Antibody Conjugate to 990 µL of 1X Diluent for each test strip to be used for testing (1/100 dilution). Dilute immediately before use and protect from light. Mix uniformly, but gently. Avoid foaming.

### 9.4 Antibody coated plate

Ready to use as supplied. Unseal foil pouch and remove plate from pouch. Remove all strips and wells that will not be used in the assay and place back in pouch and re-seal along with desiccant.

### 9.5 Calibrator (Lyophilized)

Prepare according to Section 12.

### 9.6 Chromogen Substrate and Stop Solution:

Ready to use as supplied.

## 10. Sample Collection and Storage

- All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions when handling and disposing.
- If blood samples are clotted, grossly hemolyzed, lipemic, or the integrity of the sample is of concern, make a note and interpret results with caution.
- The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.
- Known interfering substances - Azide and thimerosal at concentrations higher than 0.1% inhibits the enzyme reaction.

### 10.1 Serum:

Blood should be collected by venipuncture. The serum should be separated from the cells after clot formation by centrifugation. Remove serum and assay immediately or aliquot and store samples at  $-80^{\circ}\text{C}$  (preferably) or  $-20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

### 10.2 Plasma:

Blood should be collected into a container with an anticoagulant and then centrifuged. Assay immediately or aliquot and store samples at  $-80^{\circ}\text{C}$  (preferably) or  $-20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

### 10.3 Urine

Collect mid-stream using sterile or clean urine collector. Centrifuge to remove cell debris. Assay immediately or aliquot and store samples at  $-80^{\circ}\text{C}$  (preferably) or  $-20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

## 11. Sample Preparation

### General Sample information:

The assay requires that each test sample be diluted before use. All samples should be assayed in duplicate each time the assay is performed. The recommended dilutions are only suggestions. Dilutions should be based on the expected concentration of the unknown sample such that the diluted sample falls within the dynamic range of the standard curve. If unsure of sample level, a serial dilution with one or two representative samples before running the entire plate is highly recommended.

#### 11.1 Serum

Recommended starting dilution is 1/4. To prepare a 1/4 dilution of a sample, transfer 150  $\mu$ L of sample to 450  $\mu$ L of 1X diluent. This gives you a 1/4 dilution. Mix thoroughly.

#### 11.2 Plasma

Recommended starting dilution is 1/4. To prepare a 1/4 dilution of a sample, transfer 150  $\mu$ L of sample to 450  $\mu$ L of 1X diluent.

#### 11.3 Milk samples

Recommended starting dilution is 1/400. To prepare a 1/400 dilution of a sample, transfer 5  $\mu$ L of sample to 495  $\mu$ L of 1X diluent to achieve a 1/100 dilution. Next, dilute the 1/100 by transferring 150  $\mu$ L into 450  $\mu$ L of 1X diluent to achieve a 1/400 dilution.

## 12. Standard Curve preparation

- A standard curve should be generated each time the test is performed.

12.1 Add 1.0 ml of distilled or de-ionized water to the Calibrator and mix gently until dissolved to obtain the concentration shown on the label (the reconstituted Calibrator should be aliquoted and frozen if future use is intended) The standards need to be prepared immediately prior to use (see table below).

12.2 Mix well between each step and avoid foaming.

Note: The table below is an **example** of when the calibrator concentration is 1.520 µg/ml. Actual concentration of the standards should be calculated according to the concentration described on the vial.

Standard	Volume added to 1X Diluent (µl)	Volume of 1X Diluent	Std conc. ng/ml
6	100 µl of 1.52 µg/ml Calibrator	660 µl	200
5	300 µl std #6	300 µl	100
4	300 µl std #5	300 µl	50
3	300 µl std #4	300 µl	25
2	300 µl std #3	300 µl	12.5
1	300 µl std #2	300 µl	6.25
0	---	600 µl	0

### 13. Assay Procedure

- We recommend that you assay all standards, controls and samples in duplicate.

13.1 The Standards and the test sample(s) should be loaded into the ELISA wells as quickly as possible to avoid a shift in OD readings. Using a multichannel pipette would reduce this occurrence.

Pipette 100  $\mu$ L of the below standards in duplicate:

Standard	Concentration (ng/mL)
0	0.0
1	6.25
2	12.5
3	25
4	50
5	100
6	200

- 13.2 Pipette 100  $\mu$ L of sample (in duplicate) into pre designated wells.
- 13.3 Incubate the microtiter plate at room temperature for thirty ( $30 \pm 2$ ) minutes. Keep plate covered and level during incubation.
- 13.4 Following incubation, aspirate the contents of the wells.
- 13.5 Completely fill each well with appropriately diluted Wash Solution and aspirate. Repeat three times, for a total of four washes. If washing manually: completely fill wells with wash buffer, invert the plate then pour/shake out the contents in a waste container. Follow this by sharply striking the wells on absorbent paper to remove residual buffer. Repeat 3 times for a total of four washes.
- 13.6 Pipette 100  $\mu$ L of appropriately diluted Enzyme-Antibody Conjugate to each well. Incubate at room temperature for thirty ( $30 \pm 2$ ) minutes. Keep plate covered in the dark and level during incubation.
- 13.7 Wash and blot the wells as described in 13.4 - 13.5.
- 13.8 Pipette 100  $\mu$ L of Chromogen Substrate Solution into each well.

- 13.9 Incubate in the dark at room temperature for precisely ten (10) minutes.
- 13.10 After ten minutes, add 100  $\mu$ L of Stop Solution to each well.
- 13.11 Determine the absorbance (450 nm) of the contents of each well within 30 minutes. Calibrate the plate reader to manufacturer's specifications.

## 14. Calculations

- 14.1 Subtract the average background value (Average absorbance reading of Standard zero) from the test values for each sample.
- 14.2 Average the duplicate readings for each standard and use the results to construct a Standard Curve. Construct the standard curve by reducing the data using computer software capable of generating a four parameter logistic curve fit. A second order polynomial (quadratic) or other curve fits may also be used; however, they will be a less precise fit of the data.
- 14.3 Interpolate test sample values from standard curve. Correct for sera dilution factor to arrive at the LF concentration in original samples.

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
Printed: Tuesday, September 13, 2022, 09:27:33 Page 4 of 9

Current Workbook Path: D:\lcquan\aminoglycoside\2022\20220913-JN-genneo-KANA

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA4, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\NF-STD4, Comments: , ID: 1


Kanamycin

NF-STD4 - m/z= 163.15-163.15 RT: 1.90 - 9.90  
F: + c ESI SRM ms2 485.300 [102.179-102.181 ...


Gentamicin

NF-STD4 - m/z= 322.20-322.20 RT: 1.99 - 9.99  
F: + c ESI SRM ms2 478.325 [112.208-112.210 ...


Neomycin

NF-STD4 - m/z= 163.15-163.15 RT: 0.89 - 10.01  
F: + c ESI SRM ms2 615.325 [161.175-161.177 ...


Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA4, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\NF-STD4, Comments: , ID: 1


#2475 RT: 5.90 NL: 5.60E3  
F: + c ESI SRM ms2 485.300 [102.179-102.181 ...




#2516 RT: 5.99 NL: 1.08E3  
F: + c ESI SRM ms2 478.325 [112.208-112.210 ...

#2471 RT: 5.89 NL: 6.12E2  
F: + c ESI SRM ms2 615.325 [161.175-161.177 ...

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA4, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\NF-STD4, Comments: , ID: 1

File Name	Component Name	Calculated Conc	Units	Area	RT	Calibration Equation	Calibration Curve r^2
NF-STD4	Kanamycin	54.314	ug/kg	63696	5.90	Y = -27.5919+1173.23*X	0.9955
NF-STD4	Gentamicin	52.126	ug/kg	7978	5.99	Y = 167.119+149.856*X	0.9983
NF-STD4	Neomycin	51.951	ug/kg	5787	5.89	Y = 302.619+105.567*X	0.9956

There is no signature data to report.

**Quantify Sample Report**


Printed: Tuesday, September 13, 2022, 09:27:33 Page 6 of 9

Current Workbook Path: D:\lcquan\aminoglycoside\2022\202200913-JN-genneo-KANA

**Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA6, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\BLANK, Comments: , ID: 1**


**Kanamycin**

BLANK - m/z= 163.15-163.15 RT: 1.75 - 9.75  
F: + c ESI SRM ms2 485.300 [102.179-102.181 ...


**Gentamicin**

BLANK - m/z= 322.20-322.20 RT: 1.76 - 9.76  
F: + c ESI SRM ms2 478.325 [112.208-112.210 ...

**Neomycin**


BLANK - m/z= 163.15-163.15 RT: 0.77 - 10.01  
F: + c ESI SRM ms2 615.325 [161.175-161.177 ...

**Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA6, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\BLANK, Comments: , ID: 1**


#2412 RT: 5.75 NL: 4.04E2

F: + c ESI SRM ms2 485.300 [102.179-102.181 ...


#2418 RT: 5.76 NL: 1.57E3

F: + c ESI SRM ms2 478.325 [112.208-112.210 ...


#2422 RT: 5.77 NL: 5.28E2

F: + c ESI SRM ms2 615.325 [161.175-161.177 ...


**Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA6, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\BLANK, Comments: , ID: 1**

**Kanamycin**


**Gentamicin**

$Y = 167.119 + 149.856 * X$   $R^2 = 0.9983$  W: Equal

**Neomycin**

$Y = 302.619 + 105.567 * X$   $R^2 = 0.9956$  W: Equal

File Name	Component Name	Calculated Conc	Units	Area	RT	Calibration Equation	Calibration Curve r <sup>2</sup>
BLANK	Kanamycin	NC	ug/kg	N/A	N/A	$Y = -27.5919 + 1173.23 * X$	0.9955
BLANK	Gentamicin	NC	ug/kg	N/A	N/A	$Y = 167.119 + 149.856 * X$	0.9983
BLANK	Neomycin	NC	ug/kg	N/A	N/A	$Y = 302.619 + 105.567 * X$	0.9956

There is no signature data to report.

Quantify Sample Report


Printed: Tuesday, September 13, 2022, 09:27:33 Page 7 of 9

Current Workbook Path: D:\lcquan\aminoglycoside\2022\20220913-JN-genneo-KANA

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA7, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\JN-BENDI, Comments: , ID: 1


Kanamycin

JN-BENDI - m/z= 163.15-163.15 RT: 1.75 - 9.75  
F: + c ESI SRM ms2 485.300 [102.179-102.181 ...


Gentamicin

JN-BENDI - m/z= 322.20-322.20 RT: 1.76 - 9.76  
F: + c ESI SRM ms2 478.325 [112.208-112.210 ...


Neomycin

JN-BENDI - m/z= 163.15-163.15 RT: 0.77 - 10.01  
F: + c ESI SRM ms2 615.325 [161.175-161.177 ...


Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA7, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\JN-BENDI, Comments: , ID: 1


#2412 RT: 5.75 NL: 1.49E2  
F: + c ESI SRM ms2 485.300 [102.179-102.181 ...




#2418 RT: 5.76 NL: 5.80E2  
F: + c ESI SRM ms2 478.325 [112.208-112.210 ...

#2422 RT: 5.77 NL: 1.79E2  
F: + c ESI SRM ms2 615.325 [161.175-161.177 ...

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA7, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\JN-BENDI, Comments: , ID: 1

File Name	Component Name	Calculated Conc	Units	Area	RT	Calibration Equation	Calibration Curve r^2
JN-BENDI	Kanamycin	NC	ug/kg	N/A	N/A	Y = -27.5919+1173.23*X	0.9955
JN-BENDI	Gentamicin	NC	ug/kg	N/A	N/A	Y = 167.119+149.856*X	0.9983
JN-BENDI	Neomycin	NC	ug/kg	N/A	N/A	Y = 302.619+105.567*X	0.9956

There is no signature data to report.

Quantify Sample Report


Printed: Tuesday, September 13, 2022, 09:27:34 Page 8 of 9

Current Workbook Path: D:\lcquan\aminoglycoside\2022\20220913-JN-genneo-KANA

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA7, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\JN-BEND2, Comments: , ID: 1


Kanamycin

JN-BEND2 - m/z= 163.15-163.15 RT: 1.75 - 9.75  
F: + c ESI SRM ms2 485.300 [102.179-102.181 ...


Gentamicin

JN-BEND2 - m/z= 322.20-322.20 RT: 1.76 - 9.76  
F: + c ESI SRM ms2 478.325 [112.208-112.210 ...

Neomycin


JN-BEND2 - m/z= 163.15-163.15 RT: 1.76 - 9.76  
F: + c ESI SRM ms2 615.325 [161.175-161.177 ...

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA7, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\JN-BEND2, Comments: , ID: 1


#2412 RT: 5.75 NL: 1.31E2

F: + c ESI SRM ms2 485.300 [102.179-102.181 ...


#2418 RT: 5.76 NL: 4.57E2

F: + c ESI SRM ms2 478.325 [112.208-112.210 ...


#2422 RT: 5.77 NL: 1.88E2

F: + c ESI SRM ms2 615.325 [161.175-161.177 ...


Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA7, Name: D:\data\2022\AMINOGLYCOSIDE\09\0910\JN-BEND2, Comments: , ID: 1

Kanamycin


Gentamicin

$Y = 167.119 + 149.856 * X$   $R^2 = 0.9983$  W: Equal

Neomycin

$Y = 302.619 + 105.567 * X$   $R^2 = 0.9956$  W: Equal

File Name	Component Name	Calculated Conc	Units	Area	RT	Calibration Equation	Calibration Curve r^2
JN-BEND2	Kanamycin	NC	ug/kg	N/A	N/A	$Y = -27.5919 + 1173.23 * X$	0.9955
JN-BEND2	Gentamicin	NC	ug/kg	N/A	N/A	$Y = 167.119 + 149.856 * X$	0.9983
JN-BEND2	Neomycin	NC	ug/kg	N/A	N/A	$Y = 302.619 + 105.567 * X$	0.9956

There is no signature data to report.

Quantify Sample Report

Printed: Tuesday, September 13, 2022, 09:27:34 Page 9 of 9

Current Workbook Path: D:\lcquan\aminoglycoside\2022\20220913-JN-genneo-KANA

Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA8, Name:  
 D:\data\2022\AMINOGLYCOSIDE\09\0910\SPIKE-JN100UL, Comments: , ID: 1

Kanamycin




SPIKE-JN100UL - m/z= 163.15-163.15 RT:  
 F: + c ESI SRM ms2 485.300 [102.179-102.181 ...

Gentamicin

SPIKE-JN100UL - m/z= 322.20-322.20 RT:  
 F: + c ESI SRM ms2 478.325 [112.208-112.210 ...


Neomycin

SPIKE-JN100UL - m/z= 163.15-163.15 RT:  
 F: + c ESI SRM ms2 615.325 [161.175-161.177 ...


Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA8, Name:  
 D:\data\2022\AMINOGLYCOSIDE\09\0910\SPIKE-JN100UL, Comments: , ID: 1


#2475 RT: 5.90 NL: 9.55E3  
 F: + c ESI SRM ms2 485.300 [102.179-102.181 ...

#2530 RT: 6.03 NL: 1.49E3  
 F: + c ESI SRM ms2 478.325 [112.208-112.210 ...


#2478 RT: 5.90 NL: 6.77E2  
 F: + c ESI SRM ms2 615.325 [161.175-161.177 ...


Date: 10-Sep-2022, Instrument: TSQ Quantiva, Vial: RA8, Name:  
 D:\data\2022\AMINOGLYCOSIDE\09\0910\SPIKE-JN100UL, Comments: , ID: 1

Kanamycin


Gentamicin

$Y = 167.119 + 149.856 * X$   $R^2 = 0.9983$  W: Equal

Neomycin

$Y = 302.619 + 105.567 * X$   $R^2 = 0.9956$  W: Equal

File Name	Component Name	Calculated Conc	Units	Area	RT	Calibration Equation	Calibration Curve r^2
SPIKE-JN100UL	Kanamycin	98.491	ug/kg	115526	5.90	$Y = -27.5919 + 1173.23 * X$	0.9955
SPIKE-JN100UL	Gentamicin	109.390	ug/kg	16560	6.03	$Y = 167.119 + 149.856 * X$	0.9983
SPIKE-JN100UL	Neomycin	99.361	ug/kg	10792	5.90	$Y = 302.619 + 105.567 * X$	0.9956

There is no signature data to report.